

**091206T4PTT**

**PERIOPERATIVE THEATRE TECHNOLOGY LEVEL 6**

**HE/OS/TT/CR/02/6/A**

**PREVENT NOSOCOMIAL INFECTIONS**

**July/August 2024**



**TVET CURRICULUM DEVELOPMENT, ASSESSMENT AND CERTIFICATION  
COUNCIL (TVET CDACC)**

**OBSERVATION CHECKLIST**

**INSTRUCTIONS TO ASSESSOR**

1. *Assess the candidate as the practical progresses observing the critical areas*
2. *You are required to mark the practical as the candidate perform the tasks*
3. *You are required to take video clips at critical points*
4. *Ensure the candidate has a name tag and registration code at the back and front*

**This guide consists of FOUR (4) printed pages.**

**Assessor should check the guide to ascertain that all the pages are printed as indicated and that no questions or answers are missing**

## OBSERVATION CHECKLIST

<b>Candidates Name &amp; Registration No.</b>			
<b>Assessors Name &amp; Reg code</b>			
<b>Unit(s) of competence</b>	<b>PREVENT NOSOCOMIAL INFECTIONS</b>		
<b>Venue of assessment</b>			
<b>Date of assessment</b>			
<i>Please award marks as appropriate. Give a brief comment on your observation</i>			
<b>Items to be evaluated</b>	<b>Marks Allocated</b>	<b>Marks Awarded</b>	<b>Comments</b>
1. Performed social hand wash <i>(Award 1 mark or 0)</i>	1		
2. Donned PPEs: <ul style="list-style-type: none"> <li>• Gloves</li> <li>• Apron</li> <li>• Face mask</li> <li>• Gumboots</li> <li>• Surgical cap</li> </ul> <i>(Award 1 mark for each PPE , max 5 marks)</i>	5		
<b>Task 1: Pack the laparotomy set instruments.</b>			
3. Gathered all laparotomy surgical instruments. <i>(Award 2 marks or 0)</i>	2		
4. Candidate checked the instruments for: <ul style="list-style-type: none"> <li>• Functionality</li> <li>• Stains/rusting</li> <li>• Cleanliness</li> </ul> <i>(Award 1 mark or 0 per item correctly checked)</i>	3		
5. Followed the following steps when arranging and packing the laparotomy set:	9		

<ul style="list-style-type: none"> <li>• Lined the bottom of the tray with a green tray</li> <li>• Grouped similar instruments together</li> <li>• Ensured there is adequate space between instruments</li> <li>• Placed sharp instruments separate from blunt ones</li> <li>• Placed the heavy instruments at the bottom</li> <li>• Aligned instrument handles in the same direction within the tray</li> <li>• Double wrapped the set using green towels</li> <li>• Applied sterilizing tape and secured the set with a masking tape</li> <li>• Indicated name of the set, date, time and the person who packed</li> </ul> <p><i>(Award 1 mark for each item, max 9 marks)</i></p>			
<b>Task 2: Load the packed set in the autoclave machine.</b>			
<p>6. Followed these steps when loading the set:</p> <ul style="list-style-type: none"> <li>• Checked and removed any residual water in the autoclave machine</li> <li>• Placed the wrapped or packaged instruments in the autoclave.</li> <li>• Ensured that the machine was not overloaded</li> <li>• Placed heavier sets at the bottom</li> <li>• Closed and tightened the knobs</li> </ul> <p><i>(Award 1 mark for each item, max 5 marks)</i></p>	5		
<b>Task 3: Set the autoclave machine ready for sterilization.</b>			
<p>7. Set the autoclave machine by performing the following:</p> <ul style="list-style-type: none"> <li>• Choose the correct cycle as per user manual</li> <li>• Set the time for autoclaving</li> </ul>	2		

<i>(Award 1 mark or 0 for each correct step)</i>			
8. Disposed Waste to the appropriate bins <i>(Award 1 mark or 0)</i>	1		
9. Documented the performed procedures <i>(Award 2 marks or 0)</i>	1		
10. Performed social hand wash <i>(Award 1 mark or 0)</i>	1		
<b>Total Marks</b>	<b>30</b>		
<b>Percentage = <math>\frac{X}{30} \times 100</math></b>			
<b>The candidate was found to be:</b> Competent: <input type="checkbox"/> Not Yet Competent: <input type="checkbox"/> <i>(Please tick as appropriate)</i> <i>(The candidate is competent if she/he gets at least 50% or higher.)</i>			
<b>Feedback from Candidate:</b>			
<b>Feedback to Candidate:</b>			
<b>Candidate's Signature:</b>	<b>Date:</b>		
<b>Assessor's Signature:</b>	<b>Date:</b>		