

091206T4PTT

PERIOPERATIVE THEATRE TECHNOLOGY LEVEL 6

HE/OS/TT/CR/02/6/A

PREVENT NOSOCOMIAL INFECTIONS

July/August 2024



**TVET CURRICULUM DEVELOPMENT, ASSESSMENT AND CERTIFICATION
COUNCIL (TVET CDACC)**

PRACTICAL ASSESSMENT

TIME: 45 MINUTES

INSTRUCTIONS TO CANDIDATE:

- a) *Do not write on the question paper.*
- b) *A separate plain paper will be provided.*
- c) *Return this paper to assessor at the end of the assessment*

This paper consists of TWO (2) printed pages.

You should check the paper to ascertain that all the pages are printed as indicated and that no questions or answers are missing.

PRACTICAL TASKS

1. The autoclaving process is a reliable and widely used method for sterilizing medical instruments, laboratory equipment, and other items that require complete microbial eradication. In this assessment, you are required to complete the following tasks:

Task 1: Pack the laparotomy set instruments.

Task 2: Load the packed set in the autoclave machine.

Task 3: Set the autoclave machine ready for sterilization.

2. You are provided with the following resources for the tasks
 - i. Running water
 - ii. Liquid soap
 - iii. Color Coded Bins
 - iv. Gloves
 - v. Gumboots
 - vi. Fools cap
 - vii. Pen
 - viii. Dry hand towels
 - ix. Apron
 - x. Surgical Cap
 - xi. Assorted instruments from a Laparotomy set
 - xii. Autoclave machine
 - xiii. Sterilizing tape
 - xiv. Masking tape
 - xv. Two medium green towels

THIS IS THE LAST PRINTED PAGE.